

BeneFusion nDS

BeneFusion nDS ex

**Infusion Supervision
System**

Operator's Manual



© Copyright 2020-2023 Shenzhen Mindray Scientific Co., Ltd. All rights reserved.

Release date: March 2023


Revision: 9.0

Intellectual Property Statement

SHENZHEN MINDRAY SCIENTIFIC CO., LTD. (hereinafter called Mindray Scientific) owns the intellectual property rights to this product and this manual. This manual may refer to information protected by copyrights or patents and does not convey any license under the patent rights of Mindray Scientific, nor the rights of others.

Mindray Scientific intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray Scientific is strictly forbidden.

Release, amendment, reproduction, distribution, rental, adaptation, translation of this manual in any manner whatsoever without the written permission of Mindray Scientific is strictly forbidden.

mindray,  **MINDRAY**, and **BeneFusion** are the registered trademarks or trademarks owned by Mindray in China and other countries. Mindray Scientific is authorized by Mindray to use above registered trademarks or trademarks.

Responsibility on the Manufacturer Party

Contents of this manual are subject to changes without prior notice.

All information contained in this manual is believed to be correct. Mindray Scientific shall not be liable for errors contained herein nor for incidental or consequential damages in connection with the furnishing, performance, or use of this manual.

Mindray Scientific is responsible for the effects on safety, reliability and performance of this product, only if:

- all installation operations, expansions, changes, modifications and repairs of this product are conducted by Mindray Scientific authorized personnel;
- the electrical installation of the relevant room complies with the applicable national and local requirements;
- the product is used in accordance with the instructions for use.

Warranty

THIS WARRANTY IS EXCLUSIVE AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.

Exemptions

Mindray Scientific's obligation or liability under this warranty does not include any transportation or other charges or liability for direct, indirect or consequential damages or delay resulting from the improper use or application of the product or the use of parts or accessories not approved by Mindray Scientific or repairs by people other than Mindray Scientific authorized personnel.

This warranty shall not extend to

- Malfunction or damage caused by improper use or man-made failure.
- Malfunction or damage caused by unstable or out-of-range power input.
- Malfunction or damage caused by force majeure such as fire and earthquake.
- Malfunction or damage caused by improper operation or repair by unqualified or unauthorized service people.
- Malfunction of the instrument or part whose serial number is not legible enough.
- Others not caused by instrument or part itself.

Company Contact

Manufacturer: Shenzhen Mindray Scientific Co., Ltd.
Address: 6/F, Bldg 2, 1203 Nanhuan Avenue, Yutang Block, Guangming District, 518106 Shenzhen, P.R.China
Website: www.mindray.com
E-mail Address: service@mindray.com
Tel: +86 755 81888998
Fax: +86 755 26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)
Address: Eiffestraße 80, 20537 Hamburg, Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

Preface

Manual Purpose

This manual contains the instructions necessary to operate the product safely and in accordance with its function and intended use. Observance of this manual is a prerequisite for proper product performance and correct operation and ensures patient and operator safety.

This manual is based on the maximum configuration and therefore some contents may not apply to your product. If you have any question, please contact us.

This manual is an integral part of the product. It should always be kept close to the equipment so that it can be obtained conveniently when needed.

Intended Audience

This manual is geared for clinical professionals who are expected to have a working knowledge of medical procedures, practices and terminology as required for monitoring of critically ill patients.

Illustrations

All illustrations in this manual serve as examples only. They may not necessarily reflect the setup or data displayed on your equipment.

Conventions

- ***Italic text*** is used in this manual to quote the referenced chapters or sections.
- ***Bold text*** is used to indicate the screen texts.
- → is used to indicate operational procedures.

Contents

| | |
|--|--------------|
| 1 Safety | 1 - 1 |
| 1.1 Safety Information | 1 - 1 |
| 1.1.1 Warnings | 1 - 1 |
| 1.1.2 Cautions | 1 - 2 |
| 1.1.3 Notes | 1 - 3 |
| 1.2 Equipment Symbols | 1 - 3 |
| 2 Equipment Introduction | 2 - 1 |
| 2.1 Intended Use | 2 - 1 |
| 2.2 System Components | 2 - 1 |
| 2.3 Dock | 2 - 2 |
| 2.3.1 Front View | 2 - 2 |
| 2.3.2 Left View | 2 - 3 |
| 2.3.3 Rear View | 2 - 4 |
| 2.4 Screen Display | 2 - 5 |
| 2.4.1 On-screen Symbols | 2 - 6 |
| 2.4.2 Using the Pump Touchscreen | 2 - 6 |
| 2.4.3 Using the Pump On-Screen Keyboard | 2 - 7 |
| 3 Equipment Preparation | 3 - 1 |
| 3.1 Equipment Preparation Safety Information | 3 - 1 |
| 3.2 Environmental Requirements | 3 - 2 |
| 3.3 Installation | 3 - 2 |
| 3.3.1 Securing a Pump in the Dock | 3 - 3 |
| 3.3.2 Securing a Dock in the Medical Supply Unit | 3 - 3 |
| 3.3.3 Securing a Dock in the Medical Trolley | 3 - 4 |
| 3.4 Setting Up the Equipment | 3 - 4 |
| 3.4.1 Connecting the AC Mains | 3 - 4 |
| 3.4.2 Charging the Battery | 3 - 5 |
| 3.5 Setting Up the Dock | 3 - 5 |
| 3.6 Turning on the Dock | 3 - 6 |
| 3.7 Turning off the Dock | 3 - 7 |
| 4 Dock Setup | 4 - 1 |
| 4.1 Accessing the Dock Setup Menu | 4 - 1 |
| 4.2 The Network Setup Settings | 4 - 1 |
| 4.2.1 The Network Type Settings | 4 - 1 |
| 4.2.2 The LAN IP Settings | 4 - 2 |
| 4.2.3 The WLAN Settings | 4 - 2 |
| 4.2.4 The Bedside Hotspot Settings | 4 - 3 |
| 4.2.5 The WLAN IP Settings | 4 - 3 |

| | |
|--|--------------|
| 4.2.6 The Central Station Setup Settings | 4 - 4 |
| 4.2.7 The Device Discover Settings | 4 - 4 |
| 4.2.8 The ADT Setup Settings | 4 - 4 |
| 4.2.9 The Cert Manage Settings | 4 - 5 |
| 4.2.10 The HL7 Configuration Settings | 4 - 5 |
| 4.3 The Device Management Settings | 4 - 6 |
| 4.4 The Patient Information Settings | 4 - 6 |
| 4.5 Viewing the History Record | 4 - 6 |
| 4.6 Exporting the History Record | 4 - 6 |
| 4.7 The Language Settings | 4 - 7 |
| 4.8 The Alarm Settings | 4 - 7 |
| 4.9 The Night Mode Settings | 4 - 7 |
| 4.10 The Nurse Call Settings | 4 - 7 |
| 4.11 The Import and Export Settings | 4 - 8 |
| 4.12 Viewing the Version Information | 4 - 8 |
| 4.13 Viewing the Dock Information | 4 - 8 |
| 4.14 Exporting the Log | 4 - 8 |
| 4.15 The Batch Configuration Settings | 4 - 8 |
| 5 Alarms | 5 - 1 |
| 5.1 Alarm Safety Information | 5 - 1 |
| 5.2 Understanding the Alarms | 5 - 1 |
| 5.2.1 Alarm Priorities | 5 - 1 |
| 5.2.2 Alarm Indicators | 5 - 1 |
| 5.2.3 Alarm Handling Rules | 5 - 2 |
| 5.3 Setting the Alarm Sound | 5 - 3 |
| 5.3.1 Setting the Alarm Volume | 5 - 3 |
| 5.3.2 Setting the Alarm Sound Mode | 5 - 3 |
| 5.4 Nurse Call | 5 - 3 |
| 5.5 Alarm Solutions | 5 - 4 |
| 6 Managing Patient | 6 - 1 |
| 6.1 Admitting a Patient | 6 - 1 |
| 6.2 Exporting Patient Information | 6 - 1 |
| 6.3 Importing Patient Information | 6 - 1 |
| 7 Networked Communication | 7 - 1 |
| 7.1 Network Safety Information | 7 - 1 |
| 7.2 Connecting the Equipment to the CMS | 7 - 1 |
| 7.3 Connecting the Equipment to the Monitor | 7 - 2 |
| 7.4 Connecting the Equipment to the eGateway | 7 - 2 |
| 8 Maintenance | 8 - 1 |

| | |
|---|---------------|
| 8.1 Maintenance Safety Information | 8 - 1 |
| 8.2 Maintenance and Testing Schedule | 8 - 2 |
| 8.3 Testing Methods and Procedures | 8 - 2 |
| 8.3.1 Performing Visual Inspection | 8 - 2 |
| 8.3.2 Performing Power-on Test | 8 - 3 |
| 8.3.3 Checking the Battery | 8 - 3 |
| 8.4 Maintaining the Battery | 8 - 3 |
| 8.4.1 Battery Safety Information | 8 - 3 |
| 8.4.2 Installing the Battery | 8 - 4 |
| 8.4.3 Charging the Battery | 8 - 4 |
| 8.4.4 Conditioning the Battery | 8 - 5 |
| 8.5 Checking Version Information | 8 - 5 |
| 8.6 Checking the History Record | 8 - 6 |
| 8.7 Exporting the History Record | 8 - 6 |
| 8.8 Disposing of the Equipment | 8 - 6 |
| 9 Care and Cleaning | 9 - 1 |
| 9.1 Care and Cleaning Safety Information | 9 - 1 |
| 9.2 Cleaning the Dock and Infusion Supervision Base | 9 - 2 |
| 9.3 Disinfecting the Dock and Infusion Supervision Base | 9 - 2 |
| 9.4 Cleaning the Pole Clamp and Medical Trolley | 9 - 5 |
| 9.5 Disinfecting the Pole Clamp and Medical Trolley | 9 - 5 |
| 9.6 Sterilization | 9 - 6 |
| 9.7 Impact of Improper Cleaning | 9 - 6 |
| 10 Accessories | 10 - 1 |
| A Product Specifications | A - 1 |
| A.1 Classifications | A - 1 |
| A.2 Environmental Specifications | A - 1 |
| A.3 Power Supply Specifications | A - 2 |
| A.3.1 External Power Supply Specifications | A - 2 |
| A.3.2 Battery | A - 2 |
| A.4 Physical Specifications | A - 2 |
| A.5 Hardware Specifications | A - 3 |
| A.5.1 Displays | A - 3 |
| A.5.2 LEDs | A - 3 |
| A.5.3 Audio Indicator | A - 3 |
| A.5.4 Interface Specifications | A - 3 |
| A.5.5 Signal Output Specifications | A - 3 |
| A.6 Wireless Network | A - 4 |
| B EMC and Radio Regulatory Compliance | B - 1 |

B.1 EMCB - 1
B.2 Radio Regulatory ComplianceB - 6

C AbbreviationsC - 1

D Declaration of Conformity D - 1

1 Safety

1.1 Safety Information

WARNING

- Indicates a potential hazard or unsafe practice that, if not avoided, could result in death, serious injury or damage to product/property.
-
-

CAUTION

- Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor personal injury, product malfunction or damage to product/property.
-
-

NOTE

- Provides application tips or other useful information to ensure that you get the most out of the product.
-
-

1.1.1 Warnings

WARNING

- To avoid risk of electric shock, the equipment must only be connected to mains power with protective earth. If a protective earth conductor is not provided, operate it on battery power, if possible.
- To avoid explosion hazard, do not use the equipment in the presence of oxygen-rich atmospheres, flammable anesthetics, or other flammable agents.
- The equipment is not intended to be used within the Magnetic Resonance (MR) environment.
- Do not use the multiple portable socket outlets (MPSO) or AC mains extension cords. Ensure that the sum of the individual ground leakage currents does not exceed the allowable limits.

- **Do not open the equipment housings. All servicing and future upgrades must be carried out by trained and authorized personnel. Moreover, the servicing must be done only after the AC power supply is disconnected.**
 - **Do not place the equipment or accessories in any position that might cause it to fall on the patient.**
 - **Do not start an infusion unless the setup was verified to be correct.**
 - **To avoid inadvertent disconnection, route all cables in a way to prevent a stumbling hazard. Wrap and secure excess cabling to reduce risk of entanglement by patients or personnel.**
 - **Do not touch the patient and device connectors simultaneously. Otherwise leakage current may result in patient injury.**
 - **Ensure that more than two people are present to support and operate the equipment when it is transferred within the hospital, so as to prevent damage to the equipment and avoid tipping which can cause injury.**
 - **The communication distance between the Dock and BeneFusion nCS Infusion Supervision System, the Dock and BeneVision Central Monitoring System should be less than 50 m.**
 - **Support cascade for up to six shelf modules, and ensure that each shelf module is reliably fixed.**
 - **Equipments connected to the Dock must meet the requirements of IEC 60950. Only equipments designated by the manufacturer can be connected to the Dock. Based on patient safety, do not insert products that are not specified by the manufacturer into the Dock and its connectors.**
 - **To avoid electric shock, do not touch patient and other non-defibrillation proof equipments during defibrillation. Defibrillation will not affect the performance of the equipment.**
-
-

1.1.2 Cautions

CAUTION

- **Ensure that the equipment is supplied with continuous electric power during work. Sudden power failure may cause data loss.**
- **Electromagnetic fields may affect equipment performance. This makes it necessary for other equipment used in the vicinity of this equipment to meet EMC standards. Mobile phones, X ray and MRI equipment are all potential interference sources because of their high-intensity electromagnetic radiation.**
- **Always install or carry the equipment properly to avoid damage caused by drop, impact, strong vibration or other mechanical force. The equipment should be observed to verify normal operation after fall, otherwise it cannot be used.**

- Dry the equipment immediately in case of rain or water spray.
- Some settings are password protected and can only be changed by authorized personnel. Contact your department manager or biomedical engineering department for the passwords used at your facility.










1.1.3 Notes




















NOTE

- The software was developed in compliance with IEC62304.
- The equipment provides power-down storage. Alarms limit setting and history record are saved and will be maintained if the equipment is powered down suddenly. The storage time is equals to the equipment's service life. The alarm limit settings before power-down are reloaded when the equipment is restarted.
- This manual describes all features and options. Your equipment may not have all of them.

1.2 Equipment Symbols

Some symbols may not appear on your equipment.

| | | | |
|---|---|---|--|
|  | Refer to instruction manual/ booklet |  | Caution |
|  | Alternating current |  | Input/output |
|  | Battery |  | USB connector |
|  | Date of manufacture |  | Manufacturer |
| IP33 | Protected against solid foreign objects with a diameter no less than 2.5 mm in diameter. Protected against spraying liquid water. |  | DEFIBRILLATION-PROOF TYPE CF APPLIED PART |

| | | | |
|---|--|---|---|
|  | Atmospheric pressure limitations |  | Humidity limitations |
|  | THIS WAY UP |  | Keep dry |
|  | Fragile, handle with care |  | STACKING LIMIT BY NUMBER |
|  | Dispose of in accordance to your country's requirements |  | Authorized representative in the European Community |
|  | CE mark, comply with the requirements of the Council Directive 93/42/EEC (Medical Device Directive). |  | Serial number |
|  | Volume adjustment |  | Secondary Dock connector |
|  | Locking |  | Unlocking |
|  | Computer network |  | Stand-by |
|  | Non-ionizing electromagnetic radiation |  | Temperature limitations |
|  | No pushing | | |

2 Equipment Introduction

2.1 Intended Use

The Infusion Supervision System is in conjunction with the infusion pump and syringe pump, providing space management, power management, alarm management, information display, and communicate with pump to transmit data.

The Infusion Supervision System is intended for adults, pediatric and neonate.

The Infusion Supervision System is expected to be used in institutes or units with healthcare capabilities, such as operating rooms, emergency departments, wards, ICU. Nurses can easily centrally manage and operate pumps.

WARNING


- **This system is intended for use only by clinical professionals or under their guidance. It must only be used by persons who have received adequate training in its use. Anyone unauthorized or untrained must not perform any operation on it.**
-
-

2.2 System Components

The Infusion Supervision System consists of the controller and shelf module.

2.3 Dock


2.3.1 Front View



- (1) Display
- (2) Alarm light
When an alarm occurs, this lamp lights and flashes corresponding with the alarm priority:
 - High priority alarms: the lamp quickly flashes red.
 - Low priority alarms: the lamp lights in yellow without flashing.
- (3) Pump bay
Houses the pump.
- (4) Multifunctional connector
Provides power and data communication to individual pumps when pumps are secured in the Dock.
- (5) Connection rail lock
Secures the pump in place.
- (6) Shelf module
Holds the pump.

- (7) Controller
- (8) Controller handle
Lifts the shelf module. To prevent handle breakage, the controller handle can be used to carry only one shelf module containing maximum four pumps.


2.3.2 Left View



- (1) Battery LED
 - Green: the battery is being charged.
 - Flashing green: the Dock runs on battery power.
 - Off: no battery is installed, or no external power is connected when the equipment is off.
- (2) External power LED
 - On: when external power supply is connected.
 - Off: when external power supply is not connected.
- (3) Volume adjustment key
Turns down the alarm volume.

- (4) Volume adjustment key
Turns up the alarm volume.
- (5) Power switch
- (6) Unlocking knob
Turns the unlocking knob clockwise to the vertical position to remove the pump.
- (7) Alarm light
- (8) Infusion tubing guide
Secures the IV tubing.

2.3.3 Rear View




- (1) Secondary Dock connector
Communicates between the primary Dock and secondary Dock.

- (2) Multifunctional connector
 - RS232 connector
 - Nurse call connector
 - BeneLink connector
 - Connects the PCA controller
- (3) Mounting bracket
Secures the pole clamp in place. Pressing the button on the left side of the bracket allows you to remove the pole clamp.
- (4) Pole clamp
Secures the Dock to an approved IV pole. The pole clamp is adjustable to apply to IV poles of different dimensions.
- (5) AC power input connector
Connects the AC power cord.
- (6) Shelf module extension connector
Connects other shelf modules.
- (7) USB connector
Connects USB devices.
- (8) Network connector, a standard RJ45 connector
- (9) Product label

2.4 Screen Display

The following figure shows the screen:













- (1) Patient information area
Displays bed number and patient name.
- (2) System information area
Displays battery status and network status. For more information, see **2.4.1 On-screen Symbols**.

- (3) Infusion volume area
Displays the total infusion volume of the current patient within 24 hours.

2.4.1 On-screen Symbols


The following table lists the on-screen symbols:

| Symbol | Description | Symbol | Description |
|---|--|---|--|
|  | The battery works correctly. The solid portion represents the remaining charge. |  | The battery is being charged. |
|  | The battery has low power and needs to be charged. |  | The battery has critically low charge and needs to be charged immediately. Otherwise, the Dock will automatically shut down. |
|  | No battery is installed. |  | Battery fault, battery communication fault, or battery charging fault. Contact service personnel for help. |
|  | Wired network is connected. |  | Wired network is not connected. |
|  | Wireless network is connected. The solid part indicates network signal strength. |  | Wireless network is not connected. |

2.4.2 Using the Pump Touchscreen

You can use the pump touchscreen to select a screen element by pressing directly on the pump's screen.

To avoid misuse, the touchscreen is locked automatically if no operation is detected in the preset time. To manually lock the touchscreen, swipe the touchscreen from top down, and select **Lock**.




To unlock the touchscreen, select  on the pump touchscreen and swipe the slider as instructed.

NOTE

- **Wipe off any water on the touchscreen in case of rain or water spray.**
-

2.4.3 Using the Pump On-Screen Keyboard

The pump on-screen keyboard enables you to enter the pump and Dock information:

- Enter the information by selecting one character after another.
- Select the Backspace key  to delete single characters.
- Select the Caps Lock key  to switch uppercase letters and lowercase letters.
- Select the Enter key  to confirm the entry and close the on-screen keyboard.

This page intentionally left blank.

3 Equipment Preparation

3.1 Equipment Preparation Safety Information

WARNING

- Use only installation accessories specified by Mindray Scientific.
 - The equipment software copyright is solely owned by Mindray Scientific. No organization or individual shall resort to modifying, copying, or exchanging it or to any other infringement on it in any form or by any means without due permission.
 - Connect only approved devices to this equipment. Devices connected to the equipment must meet the requirements of the applicable IEC standards (e.g. IEC 60950 safety standards for information technology equipment and IEC 60601-1 safety standards for medical electrical equipment). The system configuration must meet the requirements of the IEC 60601-1 medical electrical systems standard. Any personnel who connect devices to the equipment's signal input/output port are responsible for providing evidence that the safety certification of the devices has been performed in accordance to the IEC 60601-1. If you have any questions, please contact Mindray Scientific.
 - If it is not evident from the equipment specifications whether a particular combination with other devices is hazardous, for example, due to summation of leakage currents, please consult the manufacturer or an expert in the field. A determination must be made that the proposed combination will not negatively affect the devices themselves or the patient's safety.
 - Ensure that the equipment is properly secured and positioned. Position change and severe shock can lead to changes to the delivery accuracy.
-
-

CAUTION

- The equipment should be installed by the authorized personnel.
 - Before use, verify whether the packages are intact. In case of any damage, do not apply it to patient.
-
-

NOTE

- **Save the packing case and packaging material as they can be used if the equipment must be reshipped.**
-

3.2 Environmental Requirements

The operating environment of the equipment must meet the requirements specified in this manual.

The environment where the equipment is used shall be reasonably free from noises, vibration, dust, corrosive, flammable and explosive substances. Moreover, to maintain good ventilation, the equipment shall be at least 2 inches (5cm) away from around the cabinet.

When the equipment is moved from one place to another, condensation may occur as a result of temperature or humidity differences. In this case, never start the system before the condensation evaporates.

CAUTION


- **Ensure that the equipment operating environment meets the specific requirements. Otherwise unexpected consequences, e.g. damage to the equipment, could result.**
-

3.3 Installation

NOTE

- **One pole clamp per shelf module must be used to ensure stability of the Dock.**
 - **Remove the infusion bags and tubings from the IV poles or pumps, and pumps from shelf module before transporting the Dock. Carry each component separately. Failure to do so can result in system becoming unbalanced. Two or more people are required to carry the Dock configured with multiple shelf modules.**
 - **According to IEC 60601-1, ensure that the carrying capacity of infusion support shall be four times more than the total weight of the Dock (including the controller, shelf module, and pumps). For example, if you are mounting a four-pump system with a total weight of 17.6 lbs (8 kg), the carrying capacity of the infusion support needs to be more than 70.5 lbs (32 kg).**
-

3.3.1 Securing a Pump in the Dock



(1) Unlocking knob

Before securing a pump into the Dock, ensure that the following requirements are met:

- The unlocking knob on the shelf module is in the horizontal position for the selected pump bay.
- The pole clamp is removed from the pump.
- The power cord is disconnected from the pump.


To secure the pump in the Dock, firmly push the pump until you hear that the clip engages the pump bay.

To unlock and remove the pump, hold the pump you wish to remove, then turn the unlocking knob clockwise to the vertical position and slide the pump out of the bay.

NOTE

-
- **Cascade Docks support maximum 24 pumps.**
 - **Only the BeneFusion n series and e series pumps can be secured in the Dock.**
-

3.3.2 Securing a Dock in the Medical Supply Unit



(1) Mounting bracket (2) Clip connector (3) Pole clamp handle

Connect the clip connector to the mounting bracket. Adjust the pole clamp handle to secure the IV pole.

NOTE

- **Use one pole clamp for each shelf module to ensure that the Dock is properly secured on the IV pole of the medical supply unit.**

3.3.3 Securing a Dock in the Medical Trolley

The pole clamp secures the Dock in the medical trolley. For detailed information on how to install the pole clamp and medical trolley, see *The Pole Clamp and Medical Trolley Installation Guides*.

3.4 Setting Up the Equipment

Observance of this manual is a prerequisite for proper product performance and correct operation. It ensures patient and operator safety.

3.4.1 Connecting the AC Mains

The equipment is powered by AC power supply. Before connecting the equipment to the AC mains, check the followings:

- The voltage and frequency ratings of the power line are the same as those indicated besides the AC power input.
- The both sides of the power cord connectors are free of liquid or other residue.
- The inside and surroundings of the AC power input connector are free of liquid or other residues.

To connect the AC power source, follow this procedure:

1. Connect the female end of the power cord to the AC power input.
2. Connect the male end of the power cord to a wall AC outlet.
3. Check that the external power supply indicator is on.

The external power supply indicator lies at the left side of the equipment. When the AC mains is not connected, the external power supply indicator is off. When AC mains is connected, the external power supply indicator is illuminated in green.

WARNING

- **Always use the accompanying power cord delivered with the equipment.**
 - **Before connecting the equipment to the AC mains, check that the voltage and frequency ratings of the power line are the same as those indicated on the equipment.**
 - **Use the battery if the integrity of the protective earth conductor or the protective earthing system in the installation is in doubt.**
 - **Do not touch the power connector with wet hand. Eliminate the liquid or any residue inside of or at the surroundings of the AC power input connector and power cord connectors.**
-
-

3.4.2 Charging the Battery

To optimize performance, a fully or nearly fully discharged battery should be charged as soon as possible. The battery is automatically charged when the equipment is connected to the AC power.

NOTE

- **The battery can only be charged by the equipment.**
 - **If the equipment is run by battery power, ensure that the battery is adequately charged.**
-

3.5 Setting Up the Dock

Before getting started, ensure that the pump and Dock are properly set up:

- The Dock is placed on a stable surface or properly mounted using the pole clamp, and the pump is secured in the Dock.
- The Dock is plugged in to a properly-grounded, AC power outlet. See **3.4.1 Connecting the AC Mains**.
- If the Dock is run on battery power, ensure that the battery is adequately charged.


WARNING

- **Before putting the system into operation, the operator must verify that the equipment, connecting cables and accessories are in correct working order and operating condition.**
 - **Only a shelf module containing maximum four pumps can be used on a flat surface. Larger shelf module configurations are more top heavy and have an increased risk of tipping, which could lead to patient or user injury.**
-
-

NOTE

- **Stay within 1 m of the Dock while setting it up and operating it, ensuring that you have a clear view of the interface.**
 - **The Dock use a mains plug as isolation means to the mains power. Do not locate the Dock in a place difficult to operate the mains plug.**
-

3.6 Turning on the Dock

Press the power switch  to turn on the Dock. The Dock automatically performs a self-test at startup. Check that the alarm tone is heard and the alarm lamp illuminates, one after the other, in red and yellow. This indicates that the visible and audible alarm indicators function correctly. The System Error alarm will be triggered when the self-test fails.

CAUTION


- **Check that visual and auditory alarm signals are presented correctly when the Dock is powered on. Do not use the Dock if you suspect it is not working properly, or if it is mechanically damaged. Contact your service personnel or us.**
-

NOTE

- **The Dock is automatically turned on when the external AC power supply is connected.**
-

3.7 Turning off the Dock

To turn off the Dock, follow this procedure:

1. Ensure that infusion has been completed.
2. Disconnect the line from the patient, and press  for at least three seconds to turn off the Dock.

NOTE

- **Turning off the Dock does not disconnect the Dock from the AC mains. To completely disconnect the power supply, unplug the power cord.**
-

This page intentionally left blank.

4 Dock Setup

Setting the Dock enables you to customize the Dock to best meet your needs. Accessing the **Dock Setup** menu is password protected.

This chapter describes the settings and functions in the **Dock Setup** menu.

CAUTION


- **The Dock Setup can only be changed by authorized personnel. Contact your department manager or biomedical engineering department for the passwords used at your facility.**
-

NOTE

- **Dock Setup menu and Batch Config menu are available when the pump and Dock are in proper connection.**
-

4.1 Accessing the Dock Setup Menu

To access the **Dock Setup** menu, follow this procedure:

1. Swipe the pump touchscreen from top down → select **Menu** → select **Dock Setup** → input the required password → select .
2. Select the desired tab.

4.2 The Network Setup Settings

4.2.1 The Network Type Settings

| Menu Item | Default Setting | Function |
|--------------|-----------------|---|
| Network Type | Auto | Selects what kind of network your equipment will use. Auto: the equipment automatically identify your network type. |

4.2.2 The LAN IP Settings

| Menu Item | | Default Setting | Function |
|------------------------|--|-----------------|--|
| DHCP Switch | | On | Selects whether inputting the IP Address, Subnet Mask, and Gateway is required. |
| IP Address | | 0.0.0.0 | |
| Subnet Mask | | 0.0.0.0 | |
| Gateway | | 0.0.0.0 | |
| Auto-obtain DNS Switch | | On | Selects whether inputting the IP address of Preferred DNS Server and Alternate DNS Server is required. |
| Preferred DNS Server | | 0.0.0.0 | |
| Alternate DNS Server | | 0.0.0.0 | |
| MAC Address | | / | / |

4.2.3 The WLAN Settings

| Menu Item | | Default Setting | Function |
|------------|-------------|-----------------|--|
| WLAN Setup | WLAN Band | Auto | Auto: automatically identifies the WLAN band. |
| | BGN Channel | All | Selects the type of B, G, and N channels. |
| | AN Channel | All | Selects the type of A and N channels. |
| SSID | | / | / |
| Password | | / | / |
| Security | | WEP OFF | Selects the security method. |

4.2.4 The Bedside Hotspot Settings

| Menu Item | | Default Setting | Function |
|------------|-------------|-----------------|---|
| WLAN Setup | WLAN Band | 5GHz | If you want to connect the device to a monitor, the WLAN band set on the device must be the same with that of the monitor when the monitor serves as a hotspot. |
| | BGN Channel | All channels | Selects the type of B, G, and N channels. |
| | AN Channel | All channels | Selects the type of A and N channels. |
| Password | | / | Views or changes the password of the shared hotspot. If you want to connect the device to a monitor, the shared hotspot password set on the device must be the same with that of the monitor. |

4.2.5 The WLAN IP Settings

| Menu Item | | Default Setting | Function |
|------------------------|--|-----------------|--|
| DHCP Switch | | On | Selects whether inputting the IP Address , Subnet Mask , and Gateway is required. |
| IP Address | | 0.0.0.0 | |
| Subnet Mask | | 0.0.0.0 | |
| Gateway | | 0.0.0.0 | |
| Auto-obtain DNS Switch | | On | Selects whether inputting the IP address of Preferred DNS Server and Alternate DNS Server is required. |
| Preferred DNS Server | | 0.0.0.0 | |
| Alternate DNS Server | | 0.0.0.0 | |
| MAC Address | | / | / |

4.2.6 The Central Station Setup Settings

| Menu Item | Default Setting | Function |
|----------------------------|-----------------|--|
| CMS Server Address | / | Inputs the IP addresses of the CMS you want to connect to. |
| Central Station IP Address | 0.0.0.0 | / |

4.2.7 The Device Discover Settings

Multicast helps device discovery between the Dock and Dock, the Dock and CMS. Devices in the same multicast group can be mutually discovered.

| Menu Item | Default Setting | Function |
|--------------------------|-----------------|----------|
| Multicast TTL | 1 | / |
| Multicast Address | 225.0.0.8 | |
| Master Server Address | / | |
| Master Server IP Address | 0.0.0.0 | |
| Connected Status | Disconnect | |

4.2.8 The ADT Setup Settings

If the Dock is connected with the Admit-Discharge-Transfer (ADT) server through the eGateway. You can load patient information from ADT server to the Dock.

| Menu Item | Default Setting | Function |
|----------------|-----------------|--|
| Server Address | / | Inputs the host name or IP address of the ADT gateway. |
| IP Address | 0.0.0.0 | |
| Port | 3502 | Inputs the port of the ADT gateway. |
| ADT Query | Off | Selects whether patient information can be loaded to the Dock from the ADT server. |

4.2.9 The Cert Manage Settings

| Menu Item | Default Setting | Function |
|-------------------|-----------------|--|
| Local Cert Manage | / | Delete: delete the selected certifications. |
| USB Cert Manage | / | Selects certifications you want to import from the USB memory, and then select Import: import the desired certifications from the USB memory. |

4.2.10 The HL7 Configuration Settings

You can send the realtime data, waveforms, and alarms from the Dock to the hospital servers via the HL7 protocol. This page also display the server connection status.

| Menu Item | Default Setting | Function | |
|---------------------|------------------|------------|--|
| Parameter HL7 Setup | Server Address | / | Inputs the name or IP address for the server receiving the realtime data and waveform. |
| | IP Address | 0.0.0.0 | |
| | Port | 0 | / |
| | Send Data | Off | |
| | Data Interval | 1h | |
| | Connected Status | Disconnect | |
| Alarm HL7 Setup | Server Address | / | Inputs the name or IP address for the server receiving the alarm data. |
| | IP Address | 0.0.0.0 | |
| | Port | 0 | / |
| | Send Alarms | Off | |
| | Connected Status | Disconnect | |
| SSL Encrypt | Off | / | |

4.3 The Device Management Settings

| Menu Item | Default Setting | Function |
|-------------|-----------------|---|
| Facility | / | Inputs the Facility, the Department, and the Device Name. |
| Department | | |
| Device Name | | |
| Device ID | / | Displays the device ID. |

4.4 The Patient Information Settings

| Menu Item | Default Setting | Function |
|------------------|-----------------|--|
| Data Source | Dock | Selects the data source of patient information. Dock: when the patient information of the Dock and pump is inconsistent, pump's patient information is synchronously updated as the Dock. |
| Patient ID | On | Selects whether the items can be displayed and edited from the Patient Management menu. |
| Visit Number | Off | |
| Patient Location | Fixed | <ul style="list-style-type: none">• Fixed: After a patient is discharged, only patient data are removed, the Bed No. and Room No. are retained.• Unfixed: After a patient is discharged, patient data, Bed No., and Room No. are all removed from the pump. |

4.5 Viewing the History Record

| Menu Item | Default Setting | Function |
|----------------|-----------------|---------------------------|
| History Record | / | Views the history record. |

4.6 Exporting the History Record

| Menu Item | Default Setting | Function |
|-----------------------|-----------------|-----------------------------|
| Export History Record | / | Exports the history record. |

4.7 The Language Settings

| Menu Item | Default Setting | Function |
|-----------|-----------------|---|
| Language | / | Sets the language. Note: This setting is effective after the Dock has been restarted. |

4.8 The Alarm Settings

| Menu Item | Default Setting | Function |
|-------------------------------|-----------------|--|
| Alarm Sound | Sound2 | Sets the alarm sound mode. |
| CMS/eGW Disconnected Alarm | Off | Sets whether the disconnection alarm will be triggered when the Dock is disconnected from the CMS or eGateway. |

4.9 The Night Mode Settings

| Menu Item | Default Setting | Function |
|--------------|-----------------|--|
| Night Mode | Off | Sets the night mode switch. The switch is turned on: The Dock enters night mode when the set Start Time is reached. The switch is turned off: The night mode is not available for the Dock. |
| Start Time | 18:00 | Sets the start time and end time of the night mode. |
| End Time | 7:00 | |
| Sound Volume | 2 | Sets the system volume during night mode. |

4.10 The Nurse Call Settings

| Menu Item | Default Setting | Function |
|--------------|-----------------|---|
| Switch | Off | Sets the nurse call Switch, Signal Type, Trigger Type, and Alarm Level. |
| Signal Type | Pulse | |
| Trigger Type | NORM. Open | |
| Alarm Level | High | |

4.11 The Import and Export Settings

| Menu Item | Default Setting | Function |
|----------------------|-----------------|--|
| Select Config File | / | Imports configuration file, drug library or brand library by following this procedure: connect the USB drive with the configuration file, drug library or brand library to the Dock's USB connector → select the file as needed → select Import . |
| Select Drug Library | | |
| Select Brand Library | | |
| Import | | |
| Export | / | Exports configuration file to the USB drive by following this procedure: connect the USB drive to the Dock's USB connector → select Export → input the file name → select Export . |

4.12 Viewing the Version Information

| Menu Item | Default Setting | Function |
|---------------------|-----------------|--|
| Version Information | / | Displays Software Version, Internal Version, Compile Time, Driver Software, Power Software, etc. |

4.13 Viewing the Dock Information


| Menu Item | Default Setting | Function |
|---------------|-----------------|--------------------------------------|
| USB Drive | NoDisk | Displays whether there is USB drive. |
| Cascade Setup | Primary Dock | Selects the Dock type. |

4.14 Exporting the Log

| Menu Item | Default Setting | Function |
|------------|-----------------|------------------|
| Export Log | / | Exports the log. |

4.15 The Batch Configuration Settings

To set the batch configuration, follow this procedure:

1. Swipe the pump touchscreen from top down → select **Menu** → select **User Maintenance** → input the required password → select  → select **Batch Config**.
2. Set **Therapy Config in Batches**: select the desired options.
3. Set **System Config in Batches**: select the desired options.
4. Select **Batch Assign** → select **Yes** to assign the configurations to all pumps in the same Dock.

This page intentionally left blank.

5 Alarms

5.1 Alarm Safety Information

WARNING

- A potential hazard can exist if different alarm presets and default settings are used for the same or similar equipment in the same care area, for example an intensive care unit or cardiac operating room.
 - The equipment in your care area may each have different alarm settings to suit different patients. Always check that the alarm settings are appropriate for your patient before start infusing.
 - Do not rely exclusively on the audible alarm system during an infusion. Adjustment of alarm volume to a low level may result in a hazard to the patient. Always ensure that the audio alarm volume level is adequate in your care environment. Always keep the patient under close surveillance.
 - Fully evaluate the risk before changing the alarm mode setting. New alarms may be failed to be detected if the operator is not familiar with the new sound.
 - The maximum alarm delay from the alarm status of the pump or shelf module to the alarm signal (sound and text displays) generated by the equipment should not exceed 5 sec.
 - The Network Disconnected alarm delay time from the equipment to the CMS is less than 14 sec, and other alarms delay time is less than 8 sec.
-
-

5.2 Understanding the Alarms

5.2.1 Alarm Priorities

By severity, the alarms are classified into the high priority alarms and low priority alarms.

5.2.2 Alarm Indicators

When an alarm occurs, the equipment indicates it visually and audibly. For more information, see the following table.

| Alarm priority | Alarm lamp color | Alarm lamp flashing frequency | Alarm sound interval | Alarm message | Alarm priority indicator | Duty Cycle |
|---------------------|------------------|-------------------------------|----------------------|---|--------------------------|------------|
| High priority alarm | Red | 2.0±0.6Hz | 5s (±2s) | White text or symbol inside a black box | !!! | 20% ~ 60% |
| Low priority alarm | Yellow | Not flashing | 20s (±2s) | White text or symbol inside a black box | ! | 100% |

NOTE

- **The tones of the alarm sound and the reminder sound are different.**
- **When multiple alarms occur simultaneously, the alarm messages are displayed circularly, and the sound and light of the higher priority alarm are given.**


5.2.3 Alarm Handling Rules

When an alarm occurs, the alarm screen is displayed to help you identify the problem.

If an alarm is triggered by the Dock, the alarm is expressed as audible tone and visual alarm message is displayed on the controller screen.

If an alarm is triggered by the pump, the alarm is expressed on the pump screen. For detailed information on alarms of the pump, see *BeneFusion n series and e series Pump Operator's Manuals*.

If the pump is secured in the Dock, the following rules are followed:



1. Alarm sound
 - ◆ The Dock conducts centralized control of the alarm sound of the pump, the alarm sound and interval time of the reminder sound are consistent with the pump.
 - ◆ The alarm sound of the Dock are silenced for two minutes by pressing  of the pump.
 - ◆ The alarm sound of the Dock is cleared after alarms of all the pumps are cleared.
 - ◆ The alarm sound of the pump is consistent with the Dock after the audio paused pump is secured in the Dock.

2. Alarm light

- ◆ When an alarm is triggered by the pump, the Dock and the pump produce alarm light, and the alarm lights of the Dock are consistent with the highest priority alarm of the pump.
- ◆ When an alarm is triggered by the pump, the alarm lights of the Dock are cleared after alarms of all the pumps are cleared.


5.3 Setting the Alarm Sound

5.3.1 Setting the Alarm Volume

Press  and  to set the alarm volume. The alarm volume can be set from 1 to 8, in which 1 is the minimum volume and 8 is the maximum volume, the default alarm volume is 4.

5.3.2 Setting the Alarm Sound Mode

To change the alarm mode, follow this procedure:

1. Swipe the pump touchscreen from top down → select **Menu** → select **Dock Setup** → input the required password → select  → select **Alarm**.
2. Set the **Alarm Sound**.

5.4 Nurse Call

The equipment provides a multi-function connector to output nurse call signal when a user-defined alarm occurs. To obtain nurse call signal, use the nurse call cable to connect the hospital's nurse call system with the equipment's multi-function connector.


Alarms are indicated on the nurse call device only when the following conditions are met:

- The nurse call system is enabled.
- A user-defined alarm occurs.

NOTE

- **Do not rely exclusively on the nurse call system for alarm notification. Remember that the most reliable alarm notification combines audible and visual alarm indications with the patient's clinical condition.**
-

To set the alarms that are sent to the nurse call system, follow this procedure:

1. Swipe the pump touchscreen from top down → select **Menu** → select **Dock Setup** → input the required password → select  → select **Nurse Call**.
2. Set the nurse call switch.
3. Select **Signal Type** to set the type of alarms that are sent to the nurse call system.

- ◆ **Pulse:** the nurse call signal is a pulse signal and each pulse lasts one second. When multiple alarms simultaneously occur, only one pulse signal is outputted. If an alarm occurs but the previous one is not cleared, a new pulse signal will also be outputted.
 - ◆ **Continuous:** the nurse call signal lasts until the alarm ends. That is to say the duration of a nurse call signal is equal to that of the alarm condition.
4. Select **Trigger Type** to set the work mode of the nurse call relay.
 5. Select **Alarm Level** to set the priority of alarms sent to the nurse call system.

5.5 Alarm Solutions

WARNING

- **When an alarm occurs, check the Dock's status and handle the alarm as soon as possible. If the alarms do not conform with the actual situation, contact your service personnel.**
-
-

| Alarm | Priority | Causes | Solutions |
|---------------------------|----------|--|--|
| Battery Depleted | High | The battery is depleted. | Connect the Dock to the external power source. |
| System Error | High | The Dock system faults, such as shelf module communication error, etc. | Stop using the Dock, and contact your service personnel. |
| Slave Controller Abnormal | Low | Communication of the primary Dock and secondary Dock is abnormal. | <ul style="list-style-type: none"> • Disconnect the connection of the secondary Dock. • Stop using the Dock, and contact your service personnel. |
| Low Battery | Low | Low battery. | Connect the Dock to the external power source. |
| Battery Error | Low | Battery fault, such as battery over heat, charging failure, etc. | Contact your service personnel. |
| LAN IP Address Conflict | Low | The Dock and other equipments in LAN have the same IP address. | Check the network settings. |

| Alarm | Priority | Causes | Solutions |
|--------------------------|-----------------|---|---|
| WLAN IP Address Conflict | Low | The Dock and other equipments in WLAN have the same IP address. | Check the network settings. |
| Network Disconnected | Low | <ul style="list-style-type: none"> • The Dock is disconnected from the CMS, the wireless network connection symbol disconnects. • The Dock is disconnected from eGateway, the wireless network connection symbol disconnects. | <ul style="list-style-type: none"> • Reconnect the Dock with the CMS, the wireless network connection symbol restores. • Reconnect the Dock with the eGateway, the wireless network connection symbol restores. • If the alarm persists, contact your service personnel. |

This page intentionally left blank.

6 Managing Patient

6.1 Admitting a Patient

The Dock admits a new patient in the following situations:

- If the Dock is connected to the ADT server through the eGateway, patient information is loaded to the Dock after Patient ID or Visit Number is input. Then patient information in the Dock automatically update.
- If the Dock is directly connected to the CMS, patient information in the Dock and CMS synchronously updates if either patient information is edited.
- If the Dock is connected to the BeneVison N series (except for the BeneVison N1) patient monitor, patient information in the Dock synchronously update when patient information in the monitor is edited, but patient information in the Dock cannot be edited.
- If a barcode reader is connected to the Dock, entering the patient management screen and scanning the patient's barcode enters the patient's information.

6.2 Exporting Patient Information

To export the information of the current patient to the USB drive, follow this procedure:

1. Connect the USB drive to the Dock's USB connector.
2. Swipe the pump touchscreen from top down → select **Menu** → select **Patient Management** → select **Export Patient Information**.
3. Select **OK**.

Exporting the patient information automatically discharge the patient.

6.3 Importing Patient Information

To import the patient information from the USB drive, follow this procedure:

1. Connect the USB drive to the Dock's USB connector.
2. Swipe the pump touchscreen from top down → select **Menu** → select **Patient Management** → select **Import Patient Information**.
3. Select **OK**.

This page intentionally left blank.

7 Networked Communication

The equipment can be connected to the BeneFusion nCS Infusion Supervision System and BeneVision Central Monitoring System (hereafter both referred to as “CMS”), patient monitors, and the eGateway.

7.1 Network Safety Information

CAUTION

- **Wireless network designing, deploying, debugging, and maintenance should be executed by the service personnel or authorized technicians.**
 - **Always set the wireless network according to local wireless regulations.**
 - **Data communication for all network functions must be performed within a closed network or within a virtually isolated network provided by a hospital. The hospital is responsible for ensuring the security of the virtually isolated network.**
 - **Keep network authentication information, for example password, safe, protecting the network from being accessed by unauthorized users.**
 - **Do not connect non-medical devices to the network.**
 - **If wireless network signal is poor, there may be a risk of CMS data loss.**
 - **RF interference may result in wireless network disconnection.**
 - **Disconnecting from the network may result in CMS data loss and function failure. Check the patient in case of network disconnection and solve the network problem as soon as possible.**
 - **Ensure that the IP address setting is correct. Changing the network settings may result in network disconnection. Contact your service personnel if you have any problems on setting the IP address.**
-

7.2 Connecting the Equipment to the CMS

The equipment can be connected to the CMS through wired or wireless network. When connected to the CMS, the system provides the following function:

- The equipment can transmit infusion information, alarm information, and equipment information (such as battery and network) to the CMS.
- Patient information can be synchronized between the equipment and the CMS.

- Patient can be admitted or discharged by the CMS, and patient information can be transmitted to this equipment.

For more information on the CMS, see *BeneFusion nCS Infusion Supervision System and BeneVision Central Monitoring System Operator's Manuals*.

To connect the equipment to the CMS, follow this procedure:

1. Set the **IP Address**, **Subnet Mask**, and **Gateway**. For more information, see **4.2 The Network Setup Settings**.
2. Connect the equipment to the CMS through any of the following methods:
 - ◆ Admit the equipment on the CMS. Refer to the *BeneFusion nCS Infusion Supervision System and BeneVision Central Monitoring System Operator's Manuals* for details of admitting an equipment.
 - ◆ Pair the equipment on the CMS. Refer to the *BeneVision Central Monitoring System Operator's Manual* for details of pairing an equipment.
 - ◆ Set the **CMS Server Address** in the **Dock Setup** menu, and the equipment automatically search and connected to the corresponding CMS. For the setting of **CMS Server Address**, see **4.2.6 The Central Station Setup Settings**.

NOTE

- **The equipment can communicate with the CMS only when it is properly connected the CMS. If the network is interrupted, you are not able to view the infusion information through the CMS.**
-

7.3 Connecting the Equipment to the Monitor

The equipment can be connected to the BeneVision N series (except for the BeneVision N1) patient monitor through the BeneLink module.

The equipment can transmit the infusion and alarm information to the patient monitor. On the patient monitor, you can view the infusion information from the **Integrated Devices** screen and infusion trends from the **InfusionView** screen. For the detailed information, see the *BeneVision N series Operator's Manual*.

7.4 Connecting the Equipment to the eGateway

You can connect the equipment to the eGateway to implement interaction between the equipment and the Hospital Information System (HIS) via the HL7 protocol. When connected to the eGateway, the system provides the following functions:

- The equipment can transmit infusion information and drug information to the eGateway.
- Patient information can be synchronized between the equipment and the eGateway.

- The prescription system assigns prescriptions to the Dock through the eGateway and master server. Then the Dock assigns prescriptions to the specified pumps.

This page intentionally left blank.

8 Maintenance

Regular maintenance is essential to ensure that the equipment functions properly. This chapter contains information on periodic testing and maintenance.

8.1 Maintenance Safety Information

WARNING

- **To avoid electric shock, stop using the equipment if you find the housing of the equipment has signs of broken. Contact the service personnel for help in that case.**
 - **Failure on the part of the responsible individual hospital or institution using this equipment to implement a recommended maintenance schedule may cause undue equipment failure and possible health hazards.**
 - **No modification of this equipment is allowed.**
 - **This equipment contains no user serviceable parts.**
 - **The safety checks or maintenance involving any disassembly of the equipment should be performed by professional service personnel. Otherwise, undue equipment failure and possible health hazards could result.**
 - **The service personnel must be properly qualified and thoroughly familiar with the operation of the equipment.**
-
-

CAUTION

- **The equipment and accessories shall not be served or maintained while in use with a patient.**
 - **If you discover a problem with the equipment, such as the product label falls off, contact your service personnel.**
-

NOTE

- **If needed, contact the manufacture for circuit diagrams, component part lists, descriptions, calibration instructions, or other information concerning the repair of the equipment.**
-

8.2 Maintenance and Testing Schedule

Follow the maintenance and testing schedule or local regulations to perform testing and maintenance. Make sure to clean and disinfect the equipment before taking any tests and maintenance.

The following table lists the maintenance and testing schedule:

| Test/Maintenance Item | | Recommended Frequency |
|-------------------------|--------------------|---|
| Safety Tests | | |
| Electrical safety tests | | <ul style="list-style-type: none"> • Once every three years, or if required. • When the power board is repaired or replaced. • When the main board is replaced. • When the equipment drops to the ground. |
| Other Tests | | |
| Visual inspection | | Every day, before first use. |
| Power-on test | | Each time the equipment is powered on. |
| Battery check | Functionality test | <ul style="list-style-type: none"> • When the battery is first installed. • When the battery is replaced. |
| | Performance test | Every three months or if the battery runtime reduces significantly. |

8.3 Testing Methods and Procedures

Except the following maintenance tasks, all other test and maintenance tasks should be performed by the qualified service personnel only.

- Regular check, including visual inspection and power-on test
- Battery check

If your equipment needs a safety test and performance test, contact the service personnel.

8.3.1 Performing Visual Inspection

Visually inspect the equipment before it is first used every day. If you find any signs of damage, remove the equipment from use and contact the service personnel.

Verify that the equipment meets the following requirements:

- Environment and power supply specifications are met.
- The equipment housing and display screen are free from cracks or other damages.

- The power cord is not damaged and the insulation is in good condition.
- Connectors, plugs, and cables are not damaged and kinked.
- Power cord and cable are securely connected with the equipment.

8.3.2 Performing Power-on Test

The equipment automatically performs a selftest at startup. Check the following items for the power-on test:

- The equipment powers on properly.
- The alarm system works properly.
- The equipment displays properly.

8.3.3 Checking the Battery

See steps 1 to 6 of **8.4.4 Conditioning the Battery** to check battery performance. The operating time of the battery reflects their performance directly. If the operating time of a battery is noticeably shorter than that stated in the specifications, the battery may reach its service life or malfunction. If the battery performance meets the requirement, fully charge the battery again for use or charge it to 40 – 60% for storage.

8.4 Maintaining the Battery

This equipment is designed to run on rechargeable Lithium-ion battery power when the external power is not available. The equipment can switch between battery power and the external power without interrupting working. If both the external power and the battery power are available, the equipment uses the external power in preference to the battery power.

8.4.1 Battery Safety Information

WARNING

- **Use only specified battery. Use of a different battery may present a risk of fire or explosion.**
- **Do not crush, drop or puncture the battery. Mechanical abuse can lead to internal damage and internal short circuits. If a battery has been dropped or banged against a hard surface, whether damage is externally visible or not, remove the battery from use and dispose of it properly.**
- **If the battery shows signs of damage or signs of leakage, replace it immediately. Use caution in removing the battery. Avoid contacting the leakage.**
- **Extremely high ambient temperature may cause battery overheat protection, resulting in equipment shutdown.**

- **The lithium-ion battery has a service life. Replace your battery when it reaches the end of its service life. Failure to replace the battery may cause serious damage to your equipment from battery overheating.**
 - **Do not open the battery, heat the battery above 60 °C, incinerate battery, or short battery terminals. They may ignite, explode, leak or heat up, causing personal injury.**
-
-

CAUTION

- **Remove the battery if it will not be used for an extended period of time.**
-

NOTE

- **Storing the battery at high temperature for an extended period of time will significantly shorten their life expectancy.**
 - **Storing the battery in a cool place can slow the aging process. Ideally the battery should be stored at 15 °C.**
-

8.4.2 Installing the Battery

The battery must only be installed by service personnel trained and authorized by Mindray Scientific. To install the battery, contact your service personnel. The battery is installed when the equipment leaves the factory.

Replace a battery in the following situations:

- The battery has visual signs of damage.
 - The battery fails.
 - The battery is aged and its runtime significantly shorter than the specification.
 - The battery service life expires.
-

CAUTION

- **Lithium batteries replaced by inadequately trained personnel could result in a hazard, such as excessive temperatures, fire or explosion.**
 - **Properly dispose of the battery according to local regulations.**
-

8.4.3 Charging the Battery

To optimize performance, a fully or nearly fully discharged battery should be charged as soon as possible. The battery is recharged automatically when the equipment is connected to AC mains power.

NOTE

- **The battery should be charged only in this equipment.**
 - **When this equipment is used with a pump, and the equipment is connected to the AC power source, the pump battery is charged automatically.**
 - **Check the battery for adequate power when the equipment runs on battery power. Charging the battery if required.**
-

8.4.4 Conditioning the Battery

The service life of a battery depends on how frequent it is used. When properly used, the lithium-ion battery has a service life of approximately three years. If improperly used, its service life can be shortened. We recommend replacing the battery every three years.

The performance of the battery deteriorates over time. You should condition the battery every two months.

To condition a battery, follow this procedure:


1. Disconnect the equipment from the patient.
2. Turn off the equipment, and connect the equipment to the external power source.
3. Allow the battery to be charged uninterrupted until it is fully charged.
4. Disconnect the equipment from the external power source, and turn on the equipment.
5. Allow the equipment to run on the battery until the battery is completely depleted and the equipment automatically shuts down.
6. Fully charge the battery again for use or charge it to 40 – 60% for storage.

NOTE

- **If the battery is not conditioned for a prolonged time, its charge indication may not be accurate and you may wrongly evaluate the remaining battery runtime.**
 - **Do not interrupt battery conditioning.**
-

8.5 Checking Version Information


To view the system Software Version, Internal Version, Compile Time, Driver Software, Power Software, and WiFi Version, follow this procedure:

1. Swipe the pump touchscreen from top down → select **Menu** → select **Dock Setup** → input the required password → select .
2. Select **Version Information**.

8.6 Checking the History Record

The **History Record** menu shows the history of the Dock activities, including the alarms and other operations.

To access the **History Record** menu, follow this procedure:


1. Swipe the pump touchscreen from top down → select **Menu** → select **Dock Setup** → input the required password → select .
2. Select **History Record**.

NOTE

- **A total loss of power has no impact on the history records stored.**
 - **Alarms are saved as events and will remain if the equipment is powered down. The time of equipment power down is also recorded as an event.**
 - **The Dock stores up to 2500 events. When the capacity is reached, earlier events will be overwritten by later ones.**
-

8.7 Exporting the History Record

To export the history record, follow this procedure:

1. Connect the USB drive to the USB connector.
2. Swipe the pump touchscreen from top down → select **Menu** → select **Dock Setup** → input the required password → select .
3. Select **Export History Record**.

8.8 Disposing of the Equipment

The service life of this equipment is ten years. Dispose of the equipment when its service life is reached. Follow local regulations regarding the disposal of such product.

WARNING

- **For disposal of parts, batteries, packaging materials, and accessories, where not otherwise specified, follow local regulations regarding disposal of hospital waste.**
-
-

9 Care and Cleaning

In this chapter we only describe cleaning and disinfection of the Dock, infusion supervision base, pole clamp, and medical trolley. For the cleaning and disinfection of other reusable accessories, refer to their instructions for use.

9.1 Care and Cleaning Safety Information

WARNING

- **Use only the approved cleaners, disinfectants and methods listed in this chapter to clean and disinfect your equipment and accessories. Warranty does not cover damage caused by unapproved substances or methods.**
 - **Do not mix disinfecting solutions, as hazardous gases may result.**
 - **We make no claims regarding the efficacy of the listed chemicals or methods as a means for controlling infection. For the method to control infection, consult your hospital's infection control officer or epidemiologist.**
 - **Be sure to turn off the system and disconnect all power cables before cleaning the equipment.**
 - **The responsible hospital or institution shall carry out all cleaning and disinfection procedures specified in this chapter.**
-
-

CAUTION

- **Turn off the equipment and remove the power cord from the equipment before cleaning and disinfecting.**
 - **Never immerse any part of the equipment or accessories in liquids or allow liquid to enter the interior of the equipment or accessories.**
 - **Any contact of cleaners or disinfectants with connectors or metal parts may cause corrosion.**
 - **Do not pour or spray any liquid directly on the equipment or accessories or permit fluid to seep into connections or openings.**
 - **If you spill liquid on the equipment or accessories, disconnect the power supply, dry the equipment, and contact your service personnel.**
 - **Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).**
-
-

- **Dilute and use the cleaners or disinfectants according to the manufacturer's instructions.**
- **Check the equipment after cleaning and disinfecting. If there is any sign of damage, remove it from use.**

9.2 Cleaning the Dock and Infusion Supervision Base

Clean the equipment on a regular basis. Before cleaning, consult your hospital's regulations.

To clean the equipment, follow this procedure:

1. Dampen a soft lint-free cloth with water or ethanol (70%).
2. Wring excess liquid from the cloth.
3. Wipe the display screen of the equipment.
4. Wipe the external surface of the equipment with the damp cloth, avoiding the connectors and metal parts.
5. Dry the surface with a clean cloth. Allow the equipment air dry in a ventilated and cool place.

CAUTION

- **Any contact of cleaners or disinfectants with connectors or metal parts may cause corrosion.**

9.3 Disinfecting the Dock and Infusion Supervision Base

Disinfect the equipment as required in your hospital's servicing schedule. Cleaning the equipment before disinfecting is recommended. Always dilute and use disinfectants according to the manufacturer's instructions. The following table lists approved disinfectants:

| Product Name | Product Type | Manufacturer |
|---|--------------|---|
| Alpet® D2 Surface Sanitizing Wipes | Wipes | BEST SANITIZERS INC™. |
| CIDEX® OPA | Liquid | Gilag GmbH International Advanced Sterilization products |
| Clorox Dispatch® Hospital Cleaner Disinfectant Towels with Bleach | Wipes | Clorox professional products company |

| Product Name | Product Type | Manufacturer |
|---|---------------------|--|
| Clorox Healthcare® Bleach Germicidal Wipes | Wipes | Clorox professional products company |
| Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes | Wipes | Clorox professional products company |
| Diversey Oxivir® TB Wipes | Wipes | Diversey Inc |
| Metrex CaviCide1™ | Liquid, spray | METERX® RESEARCH |
| Metrex CaviWipes™ | Wipes | METERX® RESEARCH |
| PDI Sani-Cloth® AF3 Germicidal Disposable Wipe | Wipes | PDI Inc. |
| PDI Sani-Cloth® Bleach Germicidal Disposable Wipe | Wipes | PDI Inc. |
| PDI Sani-Cloth® HB Germicidal Disposable Wipe | Wipes | PDI Inc. |
| PDI Sani-Cloth® Plus Germicidal Disposable Cloth | Wipes | PDI Inc. |
| PDI Super Sani-Cloth® Germicidal Disposable Wipe | Wipes | PDI Inc. |
| VIRAGUARD® Hospital Surface Disinfectant Towelette | Wipes | VERIDIEN corporation |
| Virex® II 256 (1:256) | Liquid | Diversey Inc |
| Virex® TB | Liquid, spray | Diversey Inc |
| JIAN ZHI SU Disinfectant Tablets | Tablet | Beijing ChangJiangMai Medical Science Technology Co. Ltd |
| JIAN ZHI SU Surface Disinfectant Spray | Liquid, spray | Beijing ChangJiangMai Medical Science Technology Co. Ltd |
| JIAN ZHI SU Disinfectant, Double-chain Quaternary Ammonium | Liquid | Beijing ChangJiangMai Medical Science Technology Co. Ltd |
| DIAN'ERKANG Surface Wipes | Wipes | Shanghai Likang Disinfectant Hi-Tech Co., Ltd |

| Product Name | Product Type | Manufacturer |
|--------------------------------------|---------------------|--|
| DIAN'ERKANG Surface Disinfectant | Liquid | Shanghai Likang Disinfectant Hi-Tech Co., Ltd |
| DIAN'ERKANG Disinfectant Spray | Liquid, spray | Shanghai Likang Disinfectant Hi-Tech Co., Ltd |
| Clinell® Universal Wipes | Wipes | GAMA Healthcare Ltd |
| Clinell ® Sporicidal Wipes | Wipes | GAMA Healthcare Ltd |
| Tristel Duo™ | Liquid, foam | Tristel solutions Limited |
| Tristel Jet | Liquid, spray | Tristel solutions Limited |
| Tristel Fuse For Surfaces, 196ppm | Liquid | Tristel solutions Limited |
| Surfanios Premium, 0.25% | Liquid | ANIOS LABORATORIES |
| Surfa 'safe | Liquid, spray | ANIOS LABORATORIES |
| Wip' Anios premium | Wipes | ANIOS LABORATORIES |
| Aniosurf ND premium, 0.25% | Liquid | ANIOS LABORATORIES |
| Mikrobac® Tissues | Wipes | BODE Chemie GmbH |
| Cleanisept® Wipes | Wipes | Dr. Schumacher GmbH |
| mikrozid® PAA Wipes | Wipes | Schülke & Mayr GmbH |
| mikrozid® Sensitive Wipes | Wipes | Schülke & Mayr GmbH |
| Ecolab Incidin® OxyWipe S | Wipes | Ecolab Deutschland GmbH |
| Glutaraldehyde, 2% | Liquid | / |
| Ethanol, 70% | Liquid | / |
| Isopropanol, 70% | Liquid | / |
| Sodium hypochlorite bleach, 0.5% | Liquid | / |

| Product Name | Product Type | Manufacturer |
|--|--------------|-------------------------|
| Hydrogen peroxide, 3% | Liquid | / |
| Rely+On™ Virkon® High Level surface Disinfectant, 1% | Powder | Antec International Ltd |
| 1-Propanol, 50% | Liquid | / |
| Descosept® forte | Liquid | Dr. Schumacher GmbH |
| Descosept® AF | Liquid | Dr. Schumacher GmbH |
| Dismozon® plus, 0.4% | Powder | BODE Chemie GmbH |
| mikrozid® AF Wipes | Wipes | Schülke & Mayr GmbH |
| Terralin® Liquid | Liquid | Schülke & Mayr GmbH |
| Perform® Classic Concentrate OXY, 0.5% | Powder | Schülke & Mayr GmbH |

9.4 Cleaning the Pole Clamp and Medical Trolley

Clean the pole clamp and medical trolley on a regular basis. To clean the pole clamp and medical trolley, follow this procedure:

1. Clean the pole clamp and medical trolley with a soft cloth moistened with water or ethanol (70%).
2. Wipe off the cleaner residue with a dry cloth.
3. Allow the pole clamp and medical trolley to air dry.

9.5 Disinfecting the Pole Clamp and Medical Trolley

We recommend that the pole clamp and medical trolley should be disinfected only when necessary as determined by your hospital's policy.

Cleaning the accessories before disinfecting is recommended.

| Product Name | Product Type | Manufacturer |
|-----------------------|--------------|--------------|
| Isopropanol, 70% | Liquid | / |
| Hydrogen peroxide, 3% | Liquid | / |

| Product Name | Product Type | Manufacturer |
|--|--------------|-------------------------|
| Perform® Classic Concentrate OXY, 0.5% | Powder | Schülke & Mayr GmbH |
| Dismozon® plus, 0.4% | Powder | BODE Chemie GmbH |
| Descosept® AF | Liquid | Dr. Schumacher GmbH |
| Descosept® forte | Liquid | Dr. Schumacher GmbH |
| mikrozid® AF Wipes | Wipes | Schülke & Mayr GmbH |
| Rely+On™ Virkon® High Level surface Disinfectant, 1% | Powder | Antec International Ltd |
| Terralin® Liquid | Liquid | Schülke & Mayr GmbH |

CAUTION

- **To avoid long term damage, the accessories should be disinfected only when necessary as determined by your hospital's policy.**
-

9.6 Sterilization

Sterilization is not recommended for this equipment, related products, accessories, or supplies unless otherwise indicated in the Instructions for Use that accompany the products, accessories or supplies.

9.7 Impact of Improper Cleaning

Using cleaners other than those recommended may have the following impact:

- Product discoloration
- Metal part corrosion
- Brittle and breaking wires, connectors, and equipment housing
- Reduced cable and leadwire life
- Overall system performance degradation
- Equipment malfunction or failure

10 Accessories

The accessories listed in this chapter comply with the requirements of IEC 60601-1-2 when in use with the equipment. For details about the accessories, refer to the instructions for use provided with the accessory.

WARNING

- **Use accessories specified in this chapter. Using other accessories may cause damage to the equipment or not meet the claimed specifications.**
-
-

CAUTION

- **The accessories may not meet the performance specifications if stored or used outside the specified temperature and humidity ranges. If accessory performance is degraded due to aging or environmental conditions, contact your service personnel.**
 - **Check the accessories and their packages for any sign of damage. Do not use them if any damage is detected.**
 - **Use the accessories before the expiry date if their expiry date is indicated.**
-
-

| PN | Description |
|---------------|--|
| 0020-20-12522 | Power cord, 10A, 250V, 2.5m, International |
| 009-001075-00 | Power cord, 250V, 10A, 3m, Brazil |
| 009-001791-00 | Power cord, 250V, 16A, 3m, South Africa |
| 009-002636-00 | Power cord, 10A, 1.5m, Australia standard |
| 009-007190-00 | Power cord, 3m, India |
| 009-007191-00 | Power cord, 1.8m, Switzerland |
| DA8K-10-14452 | Power cord, USA |
| DA8K-10-14453 | Power cord, UK |
| DA8K-10-14454 | Power cord, Europe |
| 009-009837-00 | Serial port adapting cable |
| 009-009838-00 | Nurse call cable |

| PN | Description |
|---|--|
| 009-009935-00 009-014088-00 009-014606-00 | Extension Dock communication cable |
| 023-001550-00 | 2D barcode reader |
| 045-004327-00 | Pole clamp |
| 045-004155-00 | Medical trolley (for 4/8/12 pump bays) |
| 045-004356-00 | Medical trolley (for 6 pump bays) |
| 115-070311-00 | ID adapter upgrade package |
| 115-074782-00 | BeneFusion tDS Infusion Supervision Base |
| 120-022009-00 | nDS 4 pump bays (2) |
| 120-022010-00 | nDS 4 pump bays (3) |
| 120-022011-00 | nDS 4 pump bays (4) |

A Product Specifications

A.1 Classifications

The equipment is classified, according to IEC60601-1:

| | |
|--|--|
| Type of protection against electrical shock | CLASS I EQUIPMENT, equipment energized from an internal electrical power source. |
| Degree of protection against electrical shock | The pump: Defibrillation-proof type CF applied part (direct cardiac application) |
| Mode of operation | Continuous |
| Degree of protection against harmful ingress of water | IP33 |
| Degree of safety of application in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide | The equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide |
| Degree of mobility | Portable |

A.2 Environmental Specifications

| Item | Temperature (°C) | Relative humidity (noncondensing) | Barometric (kPa) |
|----------------------|------------------|-----------------------------------|------------------|
| Operating conditions | 5 to 40 | 15% to 95% | 57.0 to 107.4 |
| Storage conditions | -30 to 70 | 10% to 95% | 16.0 to 107.4 |

Storage Conditions: Corrosive-free and ventilated

WARNING

- **The equipment may not meet the performance specifications if stored or used outside the specified temperature and humidity ranges. If the performance of the equipment is degraded due to aging or environmental conditions, contact your service personnel.**
-

A.3 Power Supply Specifications

A.3.1 External Power Supply Specifications

| Item | External AC Power Supply |
|-----------|-----------------------------|
| Voltage | 100 VAC to 240 VAC |
| Current | 8A to 3.4A |
| Frequency | 50/60 Hz |
| Fuse | T2AL/AC 250V (1 controller) |

A.3.2 Battery

| | |
|-----------------------|--|
| Battery Type | Lithium-ion |
| Rated battery voltage | 7.2 VDC |
| Battery capacity | 2500mAh |
| Power supply time | Discharge for at least 2 hours (1 controller and 4 shelf modules, Wi-Fi disabled, not provide power to the pump) |
| Charging time | Charging time is not longer than 6 hours. |
| Shutdown delay | At least 30 minutes after first low battery alarm (new battery, Wi-Fi disabled, not provide power to the pump) |

A.4 Physical Specifications

| Item | Maximum Weight (kg) | L x W x H (mm) |
|-------------------------------------|---------------------|-------------------|
| Dock (1 controller and 2 pump bays) | ≤ 2.9 | ≤ 270 x 173 x 245 |
| Dock (1 controller and 4 pump bays) | ≤ 3.9 | ≤ 270 x 173 x 395 |
| Dock (1 controller and 6 pump bays) | ≤ 5.4 | ≤ 270 x 173 x 550 |
| Extension 4 pump bays | ≤ 3.0 | ≤ 270 x 173 x 335 |

A.5 Hardware Specifications

A.5.1 Displays

| Type | Size (diagonal) | Resolution |
|------|-----------------|-----------------|
| OLED | ≥ 3 inches | ≥ 256x64 pixels |

A.5.2 LEDs

| | |
|--------------------|-------------------------------------|
| Alarm lamp | 2 (two color coded: yellow and red) |
| External power LED | 1 (green) |
| Battery LED | 1 (green) |

A.5.3 Audio Indicator

| | |
|---------|--|
| Speaker | Gives alarm tones (sound pressure 50 to 68 dB). Supports multi-level tone modulation. Alarm tones comply with IEC 60601-1-8. |
|---------|--|

A.5.4 Interface Specifications

| | |
|----------------------------------|---|
| Multifunctional connector | 1 |
| Secondary Dock connector | 1 |
| USB connector | 2 |
| Network connector | 1 |
| Power input connector | 1 |
| Shelf module extension connector | 1 |

A.5.5 Signal Output Specifications

| | |
|----------------------------------|--|
| Multifunctional connector | |
| Standard | Meets the requirements of IEC 60601-1 for short-circuit protection and leakage current |
| Nurse Call Signal | |
| Driving mode | Relay drive |

| | |
|------------------------|---|
| Electric specification | ≤ 60W, ≤ 2A, ≤ 36VDC, ≤ 25VAC |
| Isolation voltage | >1500VAC |
| Action mode | Normally open or normally closed (optional) |

A.6 Wireless Network

| | |
|---|---|
| Standards | IEEE 802.11a/b/g/n/ac |
| Modulation mode | BPSK,QPSK, QAM |
| Operating frequency | 2412 MHz to 2472 MHz 5180 MHz to 5825 MHz |
| Data rate | IEEE 802.11a: 6 to 54 Mbps IEEE 802.11b: 1 to 11 Mbps IEEE 802.11g: 6 to 54 Mbps IEEE 802.11n: MCS0 to MCS7 IEEE 802.11ac: MCS0 to MCS8 |
| Transfer power | < 20 dBm (CE requirement: detection mode – RMS) < 30 dBm (FCC requirement: detection mode – PEAK) |
| Operating mode | Transmitting data through the wireless access point (AP) |
| Data security | Standards: WPA/WPA2 PSK, WPA/WPA2 EAP, WPA/WPA2 CCKM EAP methods: LEAP, EAP-TTLS, EAP-TLS, EAP-FAST, PEAP-MSChapV2, PEAP-GTC, PEAP-TLS Encryption modes: TKIP and AES |
| System capacity | Number of the Docks supported by a single AP: ≤ 16 |
| Data transmission delay between the Dock and the CMS | Total data transmission delay time between the Dock and the CMS is ≤ 8s |
| Interruption number and time between the Dock and the CMS | Total interruption duration ≤ 0.01* total communication time (Test within 24 hours, 16 Docks are roaming for 30 times) |

B EMC and Radio Regulatory Compliance

B.1 EMC

The device meets the requirements of IEC 60601-1-2: 2014.

WARNING


- **Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.**
 - **Use of this device adjacent to or stacked with other device should be avoided because it could result in improper operation. If such use is necessary, this device and the other device should be observed to verify that they are operating normally.**
 - **Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the this device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this device could result.**
 - **The non-ME EQUIPMENT (e.g. ITE) that is a part of an ME SYSTEM may be disrupted by the electromagnetic interference of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the non-ME EQUIPMENT or shielding the location.**
 - **This device is intended for use in professional healthcare facility environment only. If it is used in special environment, such as magnetic resonance imaging environment, the equipment/system may be disrupted by the operation of nearby equipment.**
-

| Guidance and Declaration - Electromagnetic Emissions | | |
|---|-------------------|---|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | | |
| Emission test | Compliance | Electromagnetic environment - guidance |
| Conducted and radiated RF EMISSIONS CISPR 11 | Group 1 | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic device. |
| Conducted and radiated RF EMISSIONS CISPR 11 | Class A | The device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic distortion IEC 61000-3-2 | Class A | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage fluctuations and flicker IEC 61000-3-3 | Complies | |

NOTE

- **The device needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.**
- **Other devices may affect this device even though they meet the requirements of CISPR.**
- **The EMISSIONS characteristics of this device make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this device might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the device.**
- **If the essential performance is lost or degraded, it may be necessary to take mitigation measures, such as re-orienting or relocating the ME EQUIPMENT or ME SYSTEM or shielding the location or stopping using the infusion pump system and contact the service personnel.**

| Guidance and Declaration - Electromagnetic Immunity | | | |
|---|---|---|--|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±8 kV contact ±15kV air | ±8 kV contact ±15kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV for input/output lines (length greater than 3 m) | ±2 kV for power supply lines ±1 kV for input/output lines (length greater than 3 m) | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ±1 kV line(s) to line(s) ±2 kV line(s) to earth | ±1 kV line(s) to line(s) ±2 kV line(s) to earth | |
| Voltage dips and Voltage interruptions IEC 61000-4-11 | 0 % U_T for 0,5 cycle 0 % U_T for 1 cycle and 70 % U_T for 25/30 cycles 0 % U_T for 250/300 cycle | 0 % U_T for 0,5 cycle 0 % U_T for 1 cycle and 70 % U_T for 25/30 cycles 0 % U_T for 250/300 cycle | Mains power quality should be that of a typical commercial or hospital environment. If the user of our product requires continued operation during power mains interruptions, it is recommended that our product be powered from an uninterruptible power supply or a battery. |
| RATED power frequency magnetic fields IEC 61000-4-8 | 30 A/m 50 Hz / 60 Hz | 30 A/m 50 Hz / 60 Hz | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| Note: U_T is the A.C. mains voltage prior to application of the test level. | | | |

| Guidance and Declaration - Electromagnetic Immunity | | | |
|--|---|------------------|--|
| The device is intended for use in the specified electromagnetic environment. The customer or the user of the device should assure that it is used in such an environment as described below. | | | |
| Immunity test | IEC 60601 Test level | Compliance level | Electromagnetic environment - guidance |
| Conducted disturbances induced by RF fields IEC61000-4-6 | 3 Vrms 150 kHz to 80 MHz | 3 Vrms | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3.5}{V} \right] \sqrt{P}$ 150k to 80 MHz |
| | 6 Vrms in ISM bands ^a between 0,15 MHz and 80 MHz | 6 Vrms | |
| Radiated RF EM fields IEC61000-4-3 | 10V/m 80 MHz to 2.7 GHz | 3V/m | $d = \left[\frac{3.5}{E} \right] \sqrt{P}$ 80 MHz to 800 MHz |
| Proximity fields from RF wireless communications equipment IEC61000-4-3 | 27 V/m 380–390 MHz | 27 V/m | $d = \left[\frac{7}{E} \right] \sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^b , should be less than the compliance level in each frequency range ^c . Interference may occur in the vicinity of equipment marked with the following symbol:  |
| | 28 V/m 430–470 MHz, 800– 960 MHz, 1700–1990 MHz, 2400– 2570 MHz | 28 V/m | |
| | 9 V/m 704–787 MHz, 5100– 5800 MHz | 9 V/m | |
| <p>Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> | | | |

^a The ISM (industrial, scientific, and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The

amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

^b Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

^c Over the frequency ranges 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended Separation Distances between Portable and Mobile RF, Communications Equipment and This Equipment

The equipment is intended for use in an electromagnetic environment in which radiated RF disturbance are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communication equipment.

| Rated Maximum Output power of Transmitter Watts (W) | Separation Distance According to Frequency of Transmitter (m) | | |
|---|--|--|---|
| | 150 kHz to 80 MHz $d = \left[\frac{3.5}{f} \right] \sqrt{P}$ | 80 MHz to 800 MHz $d = \left[\frac{3.5}{E} \right] \sqrt{P}$ | 800 MHz to 2.7 GHz $d = \left[\frac{7}{E} \right] \sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

For transmitters at a maximum output power not listed above, the recommended separation distanced in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

B.2 Radio Regulatory Compliance



The radio device used in this product is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

WARNING

- **Keep a distance of at least 20cm away from the equipment when Wi-Fi function is in use.**
-
-

C Abbreviations

| Abbreviation | In Full |
|--------------|---|
| AC | Alternating Current |
| Anti-Bolus | Anti-Bolus |
| BOLUS | Bolus |
| CCU(CICU) | Cardiac Intensive Care Unit |
| CE | Conformité Européenne |
| CISPR | International Special Committee on Radio Interference |
| CPU | Central Processing Unit |
| DC | Direct Current |
| DEFS | Dose Error Reduction Systems |
| DPS | Dynamic Pressure System |
| EEC | European Economic Community |
| EMC | Electromagnetic Compatibility |
| EMI | Electromagnetic Interference |
| EtO | Ethylene oxide |
| ICU | Intensive Care Unit |
| ID | Identification |
| IEC | International Electrotechnical Commission |

| Abbreviation | In Full |
|--------------|--|
| IEEE | Institute of Electrical and Electronic Engineers |
| ISO | International Organization for Standardization |
| IV | Intravenous |
| KVO | Keep Vein Open |
| LED | Light Emitting Diode |
| Max | Maximum |
| MDD | Medical Device Directive |
| Min | Minimum |
| MRI | Magnetic Resonance Imaging |
| N/A | Not Applied |
| OR | Operating Room |
| SN | Series Number |
| USB | Universal Serial Bus |
| VTBI | Volume To Be Infused |

D Declaration of Conformity

Declaration of Conformity V1.0



Declaration of Conformity

Manufacturer: Shenzhen Mindray Scientific Co., Ltd.
6/F, Bldg 2, 1203 Nanhuan Avenue, Yutang Block, Guangming District,
518106 Shenzhen, P. R. China

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany

Product: Infusion Supervision System

Model: BeneFusion nDS, BeneFusion nDS ex

We herewith declare that the above mentioned products meet the provisions of the Council Directive 2014/53/EU concerning radio equipment. All supporting documentations are retained under the premises of the manufacturer.

Standards Applied:

| | |
|--|---|
| <input checked="" type="checkbox"/> EN 60601-1:2006/A1:2013 | <input checked="" type="checkbox"/> EN 60601-1-2:2007/AC:2010 |
| <input checked="" type="checkbox"/> EN 62311 :2008 | <input checked="" type="checkbox"/> ETSI EN 301 489-1 V2.2.0: 2017-03 |
| <input checked="" type="checkbox"/> ETSI EN 301 489-17 V3.1.1:2017-02 | <input checked="" type="checkbox"/> EN 300 328 V2.1.1:2016-11 |
| <input checked="" type="checkbox"/> ETSI EN301 893 V2.1.1: 2017-05 | |
| <input checked="" type="checkbox"/> EN60950-1:2006+A11:2009+A1:2010+A12:2011+A2:2013 | |

Start of CE-Marking: 2020-5-15

Place, Date of Issue: Shenzhen, 2020.6.1

Signature:

Name of Authorized Signatory: Ms. Bai Yanhong

Position Held in Company: Manager, Technical Regulation

This page intentionally left blank.

