
CPR Training Manikin

# BRAVDEN

## Instructions for Use

English

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Español

Français

Deutsch

Italiano

Polski

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# INNOSONIAN


CPR Training Manikin

**BRAVDEN**

English

**Instructions for Use**

**INNOSONIAN**


## Getting started

Thank you for purchasing our product.


The purpose of CPR is to supply oxygen to the brain and heart of a casualty whose heart has stopped beating normally. Applying chest compressions correctly achieves this by causing blood to flow, supplying oxygen to the brain and heart.

'BRAYDEN' is the first CPR training manikin in the world which has a blood circulation display - designed to help the student understand and remember the principles of CPR.

The product package includes the following components.

## Product name

### ① Normal circulation indication light

This illuminates to represent that blood is reaching the brain. When the depth and speed of the compressions are correct, the light will illuminate. The brightness of the light varies according to the depth of the compressions.

### ② Blood circulation indication lights

These illuminate to represent the flow of blood from the heart to the brain. The speed of the blood flow varies according to the depth and speed of the chest compressions.

### ③ Compression depth indication lights

These illuminate to represent the volume of blood being circulated by the chest compressions. This volume depends upon the depth of the compressions - deeper compressions illuminate the lights for a longer time, shallower compressions illuminate the lights for shorter time.

### ④ Power switch

This turns the manikin on and off.

### ⑤ Power status indicator lamp

When the manikin is switched on, a green light indicates that it is operating normally, a flashing red light indicates that the batteries are running low.

### ⑥ AC adaptor connection

Connect the optional AC adapter here.

### ⑦ Clicker On/Off switch

A 'clicker' which sounds with each compression over 5cm in depth can be turned on or off using this switch.

### ⑧ Battery cover


Remove this to replace the batteries.

## How to use the manikin


### Chest compression

- Turn on the manikin.
- Ensure that the power status indicator lamp is illuminated green.

- Begin compressions. The 'blood circulation' and 'compression depth' indication lights will illuminate according to the depth and speed of compressions.
- When chest compressions with excessive depth of compression (over 6cm) are performed, all of the 'compression depth' indication lights will illuminate. When chest compressions with a fast speed (over 120 times per minute) are performed, all of the 'blood circulation' indication lights will illuminate.


- When chest compressions with a depth of compression between 5cm and 6cm and a speed between 100 times and 120 times per minute are performed, the 'normal circulation' indication light will illuminate. For more details, refer to "Change the operation method".

### Artificial respiration

- Open the face shield.

- Place the mouth part of faceshield on the mouth of manikin.

- Open the airway.

- Hold the nose with your thumb and index finger, and blow air through the mouth for 1 second.

- Visually check that the manikin chest is rising.

## Change the operation method

< Basic operation >

< Specialized operation >

- When the manikin is turned on, the manikin always operates in <basic operation>, and when chest compressions with a constant depth of compression over 3cm for 5 seconds are performed, it will switch to <specialized operation>.
- The 'normal circulation' indication light will illuminate when chest compressions with a depth of compression over 5cm and a speed between 100 times and 120 times per minute in <basic operation> or chest compressions with a depth of compression between 5cm and 6cm and a speed between 100 times and 120 times per minute in <specialized operation> are performed.

## Note

- When the manikin is turned on, a one-time automatic calibration will be carried out for approximately 2 seconds. It is recommended that chest compressions on the manikin not be performed when automatic calibration is being carried out.
- If the chest compression is not detected for 2 seconds, the blood circulation indication lights, the compression depth indication lights and the normal circulation indication light will be turned off and the manikin will enter stand-by mode. When in stand-by mode, the power status indicator lamp will flash green.
- When the user carries out chest compressions again while the manikin is in stand-by mode, the chest compressions will be detected automatically and the indicators will operate again. However, if left in stand-by mode for 30 minutes, the manikin will switch itself off to conserve the batteries. To switch the manikin on again, slide the power switch to the 'on' position for more than 2 seconds and then switch it on again.
- While in stand-by mode, the power consumption is very low. However, we recommend that you switch off the manikin when not in use.
- If the battery is low, the power status indicator lamp will flash red. At this time, replace the battery with a new one.

## Change the clicker setting

< ON >


< OFF >

- Turn the 'clicker' on and off using the switch as shown. The clicker makes a "clicking" sound when the chest is compressed more than 5cm.

## Maintenance

### Replacing the artificial lung

### Battery

#### • Replacement method

#### Caution

- Keep the used batteries out of reach of children and dispose of safely.
- Do not use new batteries and old batteries together. When replacing them, replace all 6 batteries together.
- If you do not use the manikin for a long period of time, remove and store the batteries.

## Product and accessories

### Product : CPR Training Manikin "BRAYDEN"

**IM13**  
(Blood circulation display type)

**IM13-S**  
(Simple type)

The IM13-S model does not include the blood circulation display function. If you want the blood circulation display function, you need to purchase the IM13 model.

### Product specifications (IM13)

<b>Dimensions</b>	32cm x 59cm x 21cm (W x H x D)
<b>Weight</b>	3.1Kg (6.8lb)
<b>Operating temperature</b>	0°C to 40°C (32°F to 104°F)
<b>Storage temperature</b>	-20°C to 60°C (-4°F to 140°F)
<b>Humidity</b>	Relative humidity between 5% and 95%
<b>Power</b>	AC adaptor (Input AC 100-260V, Output DC 9V 1A) Battery C type 1.5V (6 batteries)

### Accessories

Face skin  
IM13-SA01

Body skin  
IM13-SA02

Head unit  
IM13-SA03

Chin  
IM13-SA04

Body  
IM13-SA05

Rib  
IM13-SA06

Face connection  
IM13-SA07

Airway valve  
IM13-SA08

Artificial lung  
IM13-SA09

Clicker  
IM13-SA10

Clicker actuator  
IM13-SA11

Spring  
IM13-SA12

Mat bag  
IM13-SA13

Instruction for use  
IM13-SA14

Faceshield  
IM13-SA15

AC adaptor  
IM13-OA01

## Quality assurance

The manikin is guaranteed against defects for 3 years from the date of purchase. This guarantee covers the mechanical parts of the manikin only, the electronic components are guaranteed for 2 years from the date of purchase. Consumable items are excluded from this guarantee.

If making a claim, you must enclose a proof of purchase and details of where the manikin was purchased. In the event that you are not able to supply a proof of purchase, or if it is illegible, free after sales service may be restricted.

Also, in any of the following cases, after-sales service may be restricted.

- Damage or defect of the manikin due to abnormal use, user's negligence or mistake
- Damage or defect of the manikin due to repair of the manikin at a place other than the company, authorized sales agency or service center or replacement of component with other company's product which is not provided by us
- Damage or defect of the manikin due to its use or storage in unsuitable environments such as in high temperature, high humidity, or where there is existence of chemical compounds or microorganisms
- Damage or defect of the manikin due to natural disasters or acts of God such as fire, earthquake, flood and lightning and other events that are beyond the scope of the company's control
- When information such as the date of purchase, customer's name, name of supplier and manufacturer's serial number has been modified intentionally
- When there is no proof of purchase

## Service center

### Innosonian, Inc.

Homepage	www.innosonian.com
Address	4F, 62, Bangbae-ro, Seocho-gu, Seoul, 06704, Republic of Korea
Tel	+82 2-533-9410
Fax	+82 70-8673-5285

### Product purchase and service inquiry

Tel	+82 70-4801-5378
E-mail	info@innosonian.com

This manikin observes the environmental requirements specified in EMC Guideline 2004/108/EC and EU Guideline RoHS 2011/65/EU for harmful substances.

### Correct disposal of the batteries in this product (Applicable in countries with separate collection systems)


This marking on battery, manual or packaging indicates that the batteries in this product should not be disposed of with other household waste at the end of their working life. Where marked, the chemical symbols Hg, Cd or Pb indicate that the battery contains mercury, cadmium or lead above the reference levels in EC Directive 2006/66. If batteries are not properly disposed of, these substances can cause harm to human health or the environment.

To protect natural resources and to promote material reuse, please separate batteries from other types of waste and recycle them through your local, free battery return system.

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