

The materials and information provided via this brochure are distributed internationally, but not all the products or services found on this brochure are necessarily available in your country or area. Please contact local office or agents to make sure of the availability of the relevant products or services.

**Mindray Medical Netherlands B.V.**

Drs.W.van Royenstraat 8, 3871 AN Hoevelaken, The Netherlands  
Voice: (31-33) 254-4911  
Fax: (31-33) 253-4280

**Mindray Medical Italy S.r.l.**

Via Leonardo DA Vinci, 158, 20090 Trezzano Sul Naviglio (MI), Italy  
Voice: (39-02) 5737-401  
Fax: (39-02) 5523-2018

**Mindray Medical España S.L.**

Avenida Manoteras 38, Edificio B, Bajo B008, 28050-Madrid Spain  
Voice: (34-91) 392-3754  
Fax: (34-91) 392-3756

**Mindray Medical Germany GmbH**

Zwischen den Bächen 4, 64625 Bensheim, Germany  
Voice: (49) 6251-17524-0  
Fax: (49) 6251-17524-20

**Mindray (UK) Limited**

3 Percy Road, St Johns Park, Huntingdon, Cambs, PE29 65Z United Kingdom  
Voice: (44) 1480-416-840  
Fax: (44) 1480-436-588

**Mindray Medical France SARL**

Europarc Creteil 1 allée des cerisiers 94000 CRETEIL, France  
Voice: (33-1) 4513-9150  
Fax: (33-1) 4513-9151

**Mindray Medical Sweden AB**

Rissneleden 136, 174 57 Sundbyberg, Sweden  
Voice: (46) 8 555 54 121  
Fax: (46) 8 555 54 101

**Mindray Medical Rus Co.Ltd**

2 Zvenigorodskaya Street, 13, bld 41, 123022 Moscow, Russia  
Voice: (7-499) 553-6036  
Fax: (7-499) 553-6039

**Mindray Medical Colombia S.A.S**

Room No.702, Torre UNIKA,Carrera 9 No.77-67, Bogota, Colombia  
Voice: (57-1) 3130892 / 3210916

**Mindray do Brasil – Comércio e Distribuição de Equipamentos Médicos Ltda.**

Rua Tavares Bastos, 329, Perdizes, São Paulo - SP, CEP 05012-020  
Voice: (55-11) 3124-8026  
Fax: (55-11) 3078-8035

**Mindray Medical Mexico S. de R.L. de C.V.**

Félix Parra # 175 Colonia, San José Insurgentes Sur Delegación Benito Juárez 03900 México, D. F.  
Voice: (52-55) 5661-9450 / (52-55) 5662-6620  
Fax: (52-55) 5662-6597

**PT.Mindray Medical Indonesia**

The East 11th Floor, Unit 08 Jl. Dr. Ide Anak Agung Gde Agung Kav. E 3.2 No.1 Jakarta 12950  
Voice: (6221) 2902-7280  
Fax: (6221) 2902-7283

**Mindray Medical Egypt**

Ground Floor -11 El Imam Metwally El Sharawi Street, Sheraton Heliopolis - El Nozha, Cairo, Egypt  
Voice: (2) 02 22691047  
Fax: (2) 02 22691046

**Mindray Medical India Pvt. Ltd. (Delhi Office)**

16th floor, Building No. 9-B, DLF Cyber City, DLF Phase- 3, Gurgaon- 122002, Haryana  
Voice: (91-124) 4632-488  
Fax: (91-124) 4632499

**Mindray Medical India Pvt. Ltd. (Mumbai office)**

B-404, City Point, Andheri Kurla Road, Andheri East, Mumbai, 400 059, India  
Voice: (91-22) 4020-0000  
Fax: (91-22) 4020-0011

**Mindray Medikal Teknoloji Istanbul Ltd. Sti.**

KOYALTI MEVKII 29 EKIM CAD.N.1 KUYUMCUKENT 2.NO.PLZ.3.KT. D.1- 9 34520 BAHÇELIEVLER/YENIBOSNA, ISTANBUL, TÜRKİYE  
Voice: (90-212) 482 08 77  
Fax: (90-212) 482 08 78

**Mindray Medikal Teknoloji Istanbul Ltd. Sti. Branch Office Ankara**


TURAN GÜNEŞ BULVARI ALEKSANDER DÜBÇEK CAD. NO:18/A ÇANKAYA/ANKARA, TÜRKİYE  
Voice: (90-312) 441 96 92  
Fax: (90-312) 441 96 93

**Mindray Medical Vietnam Co.,Ltd.**

Unit 2, 4th Floor, "B" Building, Broadway Office Park, 102 Nguyen Luong Bang Blvd., Tan Phu Ward, Dist. 7, HCMC, Vietnam  
Voice: (84) 8 5413-6970  
Fax: (84) 8 5413-6971

**Mindray Medical Thailand Limited**

1768 Thai summit Tower Room 407, 4th floor, New patchaburee Road, Kwaeng Bangkapi, Khet Huai Khwang, Bangkok 10320, Thailand  
Voice: (66) 2652-7288  
Fax: (66) 2652-7266



**BeneView TDS**

One adaptive monitoring solution whenever wherever



|  |   |
|--|---|
| <p>Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China<br/>Tel: +86 755 8188 8998 Fax: +86 755 26582680<br/>E-mail: intl-market@mindray.com www.mindray.com<br/><b>Mindray is listed on the NYSE under the symbol "MR"</b></p> | <p><b>mindray</b> are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD.<br/>©2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.<br/>P/N:ENG-BeneView TDS-210285x8P-20130929</p> <p style="text-align: center;"><b>mindray</b></p> |
|--|---|

**mindray**  
healthcare within reach



# The only one monitoring solution to reinvent your adaptive care



## TDS handle/module rack

The TDS handle combines the flexibility of handle and BeneView module rack together. The plug and play EtCO2 and PiCCO2 extension modules work with T1 also during transport.



## BeneView T1 MPM module/transport monitor

BeneView T1 combines the MPM module and transport monitor together. It's capable of plug and play with BeneView host monitors directly. T1 provides 5" touchscreen, 5 hours battery and built-in wireless network.



## TDS Docking Station

The TDS docking station provides quick plug & dock as well as I/O extension for BeneView T1 or TDS handle. The magic flexibility guarantees better adaptive care.



It's standalone modular monitoring with 19" touchscreen.  
 It's very compact transport monitoring.  
 It's wireless mobile monitoring.  
 It's plug and play MPM modular.  
 It's companion monitoring with your bedside monitor.  
 The TDS solution combines magic things to reinvent your adaptive care.  
 It's all you need for patient monitoring to adapt the changes day to day.



## BeneView companion mode

The TDS Docking Station connects with its host monitor and works at companion mode. The patient will be monitored in T1 as well as the host BeneView simultaneously. Simply quick release the TDS handle or T1 whenever patient transferring needed.



## BeneView plug-in modules


The TDS handle is compatible with BeneView sidestream, mainstream and Microstream EtCO2 modules as well as PiCCO module. The plug and play extension modules quickly improve the monitoring capability and flexibility at anywhere anytime.



## TDS external touchscreen

The 19" touchscreen connects with TDS Docking Station and provides very intuitive display and convenient navigation for BeneView T1 at the bedside.





# Seamless information workflow across the entire care cycle

## Adaptive care everywhere

The TDS can be used throughout the points of care, from the ambulance and emergency room, over the peri-operative areas, intensive care and intermediate care units to the general ward. TDS is small enough, light enough and powerful enough to follow your patient anywhere, anytime, continuously and seamlessly. Together with the BeneView host monitor, it helps streamlining your workflow, enhance patient safety, optimize your resource management and reduce costs. It allows caregivers focusing on their patients at every step of the care process, both in- and outside the hospital.



### Ambulance

Share the patient data between ambulance and hospital for earlier diagnosis and treatment



### Emergency room

Clearly monitors patient vital signs and immediately take for transport



### Patient transport

Follows the patient at every step with very compact and comprehensive monitoring including EtCO2 and PiCCO



### Peri-operative care areas

Quickly combines with bedside BeneView for powerful anesthesia monitoring workstation



### Intensive care unit


Proceeds high-acuity monitoring with companion views and closer control



### Intermediate care areas

Seamless information from the first moment in ER with built-in 120 hours trends and events memory





# Continuous patient monitoring through and through

## Adaptive care on the go

The TDS is very compact and lightweight to be easily carried and goes with the ambulatory patient and caregiver. Its auto-switchable wireless and wired network connect seamlessly with Hypervisor central stations and CMS remote viewer. The patient data is centrally managed and remote accessible across the whole IT infrastructure with electronic medical record. It helps to liberate the caregivers' hands and improve the customer satisfaction. It enables more informed care decisions as well as faster responses.



### At point of care

Synchronizes the patient data and settings with the Hypervisor central station wirelessly



### Patient ambulatory

Frees the caregiver's hands and goes with the patient everywhere for continuous monitoring around-the-clock.



### Caregiver on the go

Automatically switch between hardwired and wireless connection during the Take-Go-Plug process



### Rolling to the bedside

Queue the historic monitoring data and upload to the central station whenever connection resumed

