

BeneFusion eSP (PCA)

Patient-controlled analgesia

Data Sheet



Physical Specifications

Weight	≤ 1.6kg
Size	≤ 257x 150 x 73mm
Screen	3.5 inch touchscreen TFT color LCD, 200x400 pixels
Brightness	1-8 levels, adjustable
Display	Drug name, Bolus Dose, Lock time, Rate, Drug Concentration, Effective/Actual times, Volume, Pressure, Battery, Time, Mode, Syringe size

PCA parameters

PCA infusion:

- Loading Dose+ PCA Bolus,
- Continuous Infusion + PCA Bolus,
- Loading Dose + Continuous Infusion + PCA Bolus,
- Single PCA Bolus

PCA Bolus Dose Unit: ml, Dose (ng, ug, mg, g, mU, U, KU, EU, mmol, mol, mcal, cal, kcal, and mEq)

PCA Bolus Dose: 0.01-99.99ml; 0.001 to 999.9 (for other units)

Lock time: 1-999min

Drug Conc.: 0.001-9999.99

Bolus Limit: 0.01ml/1h to 999.9ml/1h; 0.001 to 999.9 (for other units)

Bolus Rate: 0.01-2300ml/h

Loading Rate: 0.01-2300 ml/h

Continuous Infusion Rate:0.01-2300 ml/h

Security: with password protected screen lock, driver head lock and syringe barrel clamp lock

General Specifications

Accuracy	≤ ±1.8%
Mode	PCA Mode, Rate mode, Dose Mode, Dose Time Mode, Time mode, Sequential Mode, Intermittent Mode, Loading Dose Mode, Ramp Mode, Micro-infusion Mode
Flow rate	0.01-2300ml/h
Increment	0.01ml/h(0.01-99.99ml/h);0.1ml/h(100.0-999.9ml/h); 1ml/h(1000-2300ml/h)
Preset volume (VTBI)	0.01 ml - 9999.99 ml, increment: 0.01mL
Preset time	00:00:01-99:59:59, adjustable
Accumulated volume	0.00 ml - 99999.99 ml
KVO	0.01 to 50ml/h; increment: 0.01ml/h
Purge rate	0.01-2300ml/h
Bolus rate	0.01-2300ml/h (automatic or manual)
Occlusion detection	50-1125mmHg (15 levels selectable, respectively are 50, 150, 225, 300, 375, 450, 525, 600, 675, 750, 825, 900, 975, 1050, 1125mmHg) default is 450mmHg; 4 units of pressure selectable: mmHg/kPa/bar/psi
Anti-bolus	Unexpected bolus reduced when the occlusion occurs
Screen lock	Support password-protected screen lock
Auto-lock time	1 - 5 minutes selectable; ON/OFF switchable

Drug library	Up to 5000 drugs; 50 categories (optional); 30 care areas/departments, support color-coding drug name
History log	Up to 3500 records
Volume collection	Available in 4 methods: 24H total, current total, period, timing volume
DERS (Dose Error Reduction System)	Available; definition of dose limits; automatic alarms when reaching dose limits

Syringes

Compatibility	1/2/3/5/10/20/30/50/60ml; automatic recognition of syringe size
---------------	--

Alarms

Type	Audible and visual alarm
Level	High: Occlusion/ Syringe empty/ Syringe Disengaged/ Plunger Grippers Error/ PCA Cable Detached/ Battery Depleted/ VTBI Complete/ KVO Finish/ Relay Invalid/ System Error/ No Syringe Low: KVO Running/ Battery in Use/ Battery Error/CMS/eGW Disconnected/ Standby Time Expired/ Dock Connection Interrupted/ System Time Error/ Relay Invalid Soon/ Time Near End/ Reminder/ Low Battery/ Syringe Near Empty
Sound volume	1-8 levels selectable, default is level 4
Reminder	1-5 minutes selectable; ON/OFF switchable

Connectivity

Communication	Wireless
USB	Support drug library import, patient information loading, patient data export, history record export, calibration data import/ export
Multifunctional connector	RS232; nurse call connector; DC adapter
Integration	Synergy with Mindray patient monitor Connect with BeneVision Central Monitoring System

Battery

Operating time	≥ 5 hours at 5ml/h (≥ 11 hours at 5ml/h for smart battery)
Charging time	≤ 5 hours to full capacity (≤ 6 hours for smart battery)

Power Supply

	Voltage 100-240 V; frequency 50/60Hz; current 0.5-0.21A
--	---

Work Environment

Temperature	5-40°C for operating; -30-70 °C for storage
Relative humidity	15-95% for operating; 10-95% for storage
Atmosphere pressure	57.0-107.4 kPa for operating; 16.0-107.4 kPa for storage
Classification	Type CF, Protection Class I, IP33

www.mindray.com

P/N: ENG-BeneFusion eSP (PCA) Datasheet-210285x2P-20251230

©2019 Shenzhen Mindray Scientific Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.

mindray
healthcare within reach